

Universal Refrigerated Centrifuge



5930

3 types of biohazard-proof rotor available

High centrifugal force swinging bucket rotor compatible [RS-3020M: 3,430-3,200 xg] Conforms to International Safety Standard IEC61010-2-020

7 rotor types available to fit your needs

- 4 types of buckets available for the RS-3020M swinging bucket rotor
- High centrifugal force RS-3020M rotor (3,420 xg) usable with sealing caps, for centrifugal separation of acid-fast bacteria
- PT-745MS plate rotor (up to 3,100 xg) with shield covers that prevent contamination

Sample balance can be measured by eye (for rotors used at 5,000rpm or less)

Streamline routine work with convenient functions for easy operation

- Memory flashing function for easy spin-down (patented) automatically memorizes the duration of spin-down operations ranging from 1 to 90 seconds. The last flashing operation (1-90sec.) is automatically memorized and easily recalled simply by pressing the "Flashing" key.
- 5 available memory settings, with one-touch save and recall operation.
- Operation according to centrifugal force (g) as well as stepped centrifuging and integral centrifuging (g·sec) settings.
- Timer can be set up to 990 min/sec. When on hold, the elapsed time is displayed.
- Optional foot switch enables hands-free lid opening.

Conforms to International Safety Standard IEC61010-2-020 (except when RA-410M3 rotor is used)

- In the event of breakage, all fragments will be safely retained inside the centrifuge unit.
- Lid lock mechanism includes interlock and twin-lock hook system (right and left).
- Interlocking parts are protected against mechanical impact and employ a two-stage electrical protection circuit.
- Various sensors immediately detect the onset of any trouble during operation.


Tachometer port complies with GMP

• Easy verification of spinning state.

Sample temperature is maintained at 4°C at maximum speed of rotation of all rotors (at 25°C ambient temperature)




Variable acceleration and deceleration function Acceleration and deceleration curves are adjustable to best fit your samples and goals. Initial acceleration and deceleration speeds can be adjusted to accommodate easily perturbed samples. Density gradient centrifugation is also enabled.



Acceleration/deceleration characteristics

- Super-slow acceleration or deceleration is available for centrifugation of especially sensitive samples, to avoid resuspension.
- Suitable for centrifugation of easily resuspended samples, density gradient methods and for blood cell cleaning.
 The returned samples the acceleration/deceleration profile changes.
- The rotational speed N when the acceleration/deceleration profile changes can be set, enabling super-slow acceleration and deceleration for especially sensitive samples.




Integrator function $g \cdot sec$ value Identical applications of total centrifugal force are easily reproducible using the $g \cdot sec$ function


• The integrator function enables identical centrifugal effects to be applied. (A) = (B)



RS-3012M swinging bucket rotor As shown, the rotor is fitted with optional shield covers (No.027-9210). A secure clip system lets you easily attach and remove the shield covers.



Code No.053-5180 Maximum capacity is 750mLx4 bottles. Various tube racks can be used.



Code No.053-5220 8 micro-plates and 4 deep-well plates can be set. This bucket has no center bar so plates can be easily removed.



Code No.053-5170 8 micro-plates and 8 tissue culture plates can be set.



RS-3020M swinging bucket rotor As shown, the rotor is fitted with optional sealing caps (No.027-9300).

Rotor, Bucket, Tube Rack Table for Model 5930

Rotor Name		Nominal Capacity	Tube size (mm) diam. x length	Maximum Speed rpm	Maximum RCF xg	Bucket Code No.	Tube Rack Adapt Code No.
winging	bu	ucket rotor					
		140 x 6mL/RIA tube	10~12.3 x 65~100 (95)	3,300	2,150 *1	053-5310	055-5830
		80 x 5-7mL blood collecting tube *2	12~13.3 x 63~110 (105)	3,300	2,280	053-5310	055-5000
		100 x 10mL	12~14.2 x 76~110 (105)	3,300	2,280	053-5310	055-5840
		80(64) x 5-10mL blood collecting tube/15mL * 2	12~17.2 x 72~110 (105)	3,300	2,280	053-5310	055-0900
		80(64) x 15mL	15~17.2 x 76~110 (105)	3,300	2,280	053-5310	055-5850
		16 x 15/50mL	15mL 15~17.2 x 75~135 (125)	3,300	2.200	053-5310	055 5960
	A	16 X 15/50ML	50mL 27~36 x 79~115 (110)	3,300	2,280		055-5860
	*3	48 x 15mL Conical tube	17 x 121 (121)	3,300	2,280	053-5310	055-4660
		24 x 50mL Conical tube	30 x 117 (117)	3,300	2,280	053-5310	055-4610
		8 x 100mL	38~46.5 x 75~120 (116)	3,300	2,280	053-5310	055-5870
		4 x 250mL *4	60~61.8 x 100~141 (136)	3,300	2,280	053-5310	055-6880
		4 x 350mL	65~71 x 100~130 (122)	3,300	2,280	053-5310	055-6890
		12x/16x/20x Self-standing sputum tube *5	Cap diam. 42/ 39/ 31	3,300	*5	053-5310	*5
RS-3012M		72 x 15mL	15~17.2 x 75~110	3,300	2,280	053-5180	055-0210
		16 x 50mL	27~36.8 x 79~114	3,300	2,280	053-5180	055-0240
		48 x 15mL Conical tube	17 x 121	3,300	2,280	053-5180	055-0220
		16 x 50mL Conical tube	30 x 117	3,300	2,280	053-5180	055-0230
	6	4 x 250mL *4	60~61.8 x 100~141	3,300	2,280	053-5180	055-4490
		4 x 350mL	65~71 x 100~130	3,300	2,280	053-5180	055-5880
		4 x 500mL Conical tube *6	96 x 147	3,300	2,280	053-5180	431124
		4 x 650mL * 13	98 x 136	3,300	2,280	053-5180	_
		4 x 650mL * 14	84~86 x 100~125	3,300	2,280	053-5180	055-5890
		4 x 750mL	95~97.5 x 100~135	3,300	2,280	053-5180	—
	O	8 x Micro-plates *7	88 (M) x 124 (D) x 50 (U)	3,300	2.020 *1	053-5220	
	9	4 x Deep-well plates *7	88 (W) x 134 (D) x 50 (H)		2,020 *1		_
		8 x Micro-plates *7	02 (M) x 122 (D) x 26 (U)	0.000	770 *1	053-5170	055-5900
	O	8 x Tissue culture plates *7	92 (W) x 133 (D) x 36 (H)	2,000	770 * 1		—
RS-3020M *8		76 x 15mL	15~17.2 x 83~116	4,000	3,200 * 1	Equipped with rotor	055-0340
		48 x 15mL Conical tube	17 x 121	4,000	3,400 * 1	Equipped with rotor	055-0300
		16 x 50mL	27~35.8 x 93~117	4,000	3,310 *1	Equipped with rotor	055-0330
	0	20 x 50mL Conical tube	30 x 117	4,000	3,420	Equipped with rotor	055-0310
		4 x 500mL Conical tube *6	96 x 147	4,000	3,430	Equipped with rotor	431124
		4 x 650mL * ¹³	98 x 136	4,000	3,430	Equipped with rotor	—
		4 x 750mL	95~97.5 x 100~135	4,000	3,430	Equipped with rotor	_

Plate rotor

PT-21M	6 x Micro-plates	*7	86 (W) x 130 (D) x 62 (H)	5,000	3,100	Equipped with rotor		
P1-211VI	2 x Deep-well plates	*7	88 (W) X 130 (D) X 82 (H)				_	
PT-89M	8 x Micro-plates	*7	86 (M) × 120 (D) × 80 (U)	2,400	1,080	Equipped with rotor	055-6370	*9
	4 x Deep-well plates	*7	86 (W) x 130 (D) x 89 (H)					
PT-745MS	12 x Micro-plates	*7	86 (M) × 128 (D) × 00 (U)	4 200	3.100	Equipped with reter	055 6280	*9
	4 x Deep-well plates *7		86 (W) x 128 (D) x 90 (H)	128 (D) x 90 (H) 4,200		Equipped with rotor	055-6380	9

Fixed Angle rotor AT-508C angle rotor is a tool-less rotor.

			Tube size (mm)	AT-508C		RA-410M3		Bucket	Adaptor
Nominal Capacity		diam. x length	Maximum Speed rpm	Maximum RCF xg	Maximum Speed rpm	Maximum RCF xg	Code No.	Code No.	
8 x 12mL (plastic) *10		15.9~16.5 x 97~120	6,000	5,110	8,500	9,370	—	055-1120	
8	8 x 14mL *11		17 x 100	6,000	4,990	8,500	9,130	—	055-0550
8 x 15mL	8 x 15mL Conical tubes		17 x 121	6,000	5,270	8,500	9,530	—	055-1280
8 x 50r	8 x 50mL (plastic) *12		26.4~28.8 x 98~115	6,000	4,990	8,500	9,290	—	055-1110
8 x 50mL Conical tubes			30 x 117	6,000	5,350	8,500	9,690	—	—
Sartorius	Vivaspin 2		16.7 x 129	6,000	5,270	8,500	9,530	—	055-1280
	Vivaspin 6		16.7 x 129	6,000	5,270	8,500	9,530	—	055-1280
	Vivaspin 20		30 x 117	6,000	5,350	8,500	9,690	—	—
Millio and	Amicon Ultra-4		17 x 129	6,000	5,270	8,500	9,530	—	055-1280
Millipore	Amicon Ultra-15		30 x 121	6,000	5,350	8,500	9,690	—	—
 *1 For centrifugal force indications, please refer to the instruction manual. *2 Blood collecting tubes made by SEXISUI (Insepack II), Termo (Venoject II), Nipro (Neotube) or Becknon Dickinson (Vacutenar) can be used. *3 Figures in parentheses indicate speeds when using optional shield covers (No.027-9210). *4 175mL, 200mL and 225mL conical tubes can be used with 175 & 225mL conical tubes made by BD Falcon (No.352076). No.37585 cushions must be used with Nunc 200mL bottles (No.376813). *6 Corning No.431123 tubes require No.431124 cushions. *7 Verify plate strength before centrifuging. Plates in lower position may fracture. *8 Optional 027-9300 sealing caps can be used, except with 500mL conical tubes. *9 Facilitates removing plates from buckets. *10 Nunc plastic tubes, No.352006, 352018 and 352059. *12 Nunc plastic tubes, No.3110-0500, No.3110-0500, No.3112-0050, No.3114-0050, No.3114									

No.372655 cushions must be used with Nunc 200mL bottles (No.376813).
 * 5 Various sputum tube capacities can be used. Contact our distributor for details.

* 13 Herolab bottles, 650mL (PPCO, PC). * 14 Stainless bottles are available.

* Autoclave parameters should not be exceeded. Refer to operation and service manual.



PT-21M Plate rotor Can do ethanol precipitation at 3,100 xg.



PT-89M Plate rotor Accepts 8 micro-plates or 4 deepwell plates.



PT-745MS Plate rotor 12 micro-plates and 4 deep-well plates can be set. Ethanol precipitation can be performed at 3,100 xg.



AT-508C Fixed-angle rotor 50/15mL x 8 conical tubes This economical, multipurpose fixed-angle rotor accepts 12mL or 50mL tubes.



RA-410M3 Fixed-angle rotor 50/15mL x 8 conical tubes This multipurpose fixed-angle rotor accepts12mL or 50mL tubes.



Model 5	5930
Maximum Speed •	8,500 rpm
Maximum RCF •	9,690 xg
Maximum capacity	3,000 mL
Control system	Microprocessor control (brushless motor)
	Speed, RCF(g), Time, Temperature, Acceleration and Deceleration, 5 memory banks
Alarm display •	Lid open, Imbalance, Over speed, Rotor setting condition, Abnormally high temperature,
	Detected electrical abnormality in motor, inverter, lid interlock, speed sensor and temperature sensor
	Three selectable rates: Rapid, Slow and Super Slow
	From 200 to 8,500 rpm, 10 or 100 rpm increments
	Digital display, from 0 to 8,900 rpm, 10 rpm increments
	From 10 to 9,690 xg, 10 or 100 xg increments
	Digital display, from 0 to 10,620 xg, 10 g increments
	Digital display, from 1.00-9.99 x10° g · sec
	speed, RCF(g), time, temperature, acceleration and deceleration, 1 memory bank
Timer setting and indication •	Digital display, with HOLD function (continuous operation) and electric buzzer
	1 second to 100 seconds, 1 sec increments
	100 seconds to 990 seconds, 10 sec increments
	1 minute to 100 minutes, 1 min increments
Tomporature patting and indication	100 minutes to 990 minutes, 10 min increments digital display, -20°C to +40°C, 1°C increments
	R134a (CFC-free)
0	Single phase 50/60Hz, 110V±10% 15A, 115V±10% 15A
Fower requirements	Single phase 50/60Hz, 220V±10% 10A, 230V±10% 10A
Rated voltage and rated current	110V 11A, 115V 10A, 220V 7.0A, 230V 6.6A
0	110/115V 990W 3.6MJ/h, 220/230V 1.1kW 4.0MJ/h
	50cm x 65cm x 85cm (WxDxH) ^(*1) , 117kg
	IEC 61010-2-020 (Except when RA-410M3 rotor is used.)
i	CE Mark (AC220V and 230V versions)

(*1) Not including projections.

All specifications and external appearance are subject to change without notice.

* Use only rotors, adaptors and tube racks made by KUBOTA. Do not attempt to use rotors, adaptors or tube racks made by other manufacturers. * Follow all instructions in the relevant Instruction & Service Manual when operating centrifuges.

Kubota has acquired ISO 9001 and ISO 13485 certification.

Saf Pre	cautions the	nstruction manual carefully before you start its		Parts may be adversely affected by heat, humidity and specimens or deteriorated by aging resulting in functional failures or may even lose their safety characteristic to cause accidents.			
DISCONTINUE	Always discontine use when phenomena listed at right appear.	•A burnt smell comes out of the equipment.		IMPORTANT	To prevent occurrences of functional failures or accidents, turn off the power switch, pull out the power cable plug from the plug socket and be sure to contact your nearest Kubota Dealer, or an office of Kubota Corporation when there is no Kubota Dealer in your neighbourhood.		

KUBOTA CORPORATION

http://www.centrifuge.jp

29-9 Hongo 3-chome, Bunkyo-ku, Tokyo 113-0033, Japan Tel +81 3 3815 1331 Fax +81 3 3814 2574 E-mail: f-trade@kubotacorp.co.jp

Note: This catalogue is not for distribution in the USA, Canada, Mexico and Australia as products shown are not for sale in these countries.